EN CELLBOSTER® HAIR
Hyaluronic Acid 18 mg Strengthening Complex

INSTRUCTIONS FOR USE

DESCRIPTION

CELLBOSTER® HAIR: Hyaluronic Acid Strengthening Complex is produced using patented CHAC Technology. It consists of non-crosslinked Hyaluronic Acid with Amino Acids and Vitamins. The device is designed for strengthening hair follicles, improving hair growth and microcirculation. It restores scalp skin structure and provides hair fiber restructuring.

The device is a clear, colorless to light yellow, sterile, injectable solution for single use presented in a case of 6 x 3 mL vials.

COMPOSITION

Each vial contains:

- non-crosslinked Sodium Hyaluronate 18 mg;**
- amino acids** (Arginine, Cysteine, Glutamine, Glycine, Lysine);
- vitamins** (Niacin, Pantothenic Acid, Pyridoxine, Biotin, Cyanocobalamin, Rutin);
- Copper Gluconate, Zinc Carbonate;
- phosphate buffer pH 6.8-7.6 q.s.

INDICATIONS

The device is a resorbable injectable implant for preventing and treating various scalp and hair problems linked to hormonal imbalance, psoriasis, seborrhea, stress, nutrition and trace elements deficiencies, such as: hair loss, hair breakage, scalp dehydration, dandruff, premature graying. Improves microcirculation, restores scalp skin structure and other qualities. Provides hair fiber restructuring after physical and/or chemical damage, such as: hair dyeing, permanent wave, sun, dryer. Improves results of hair transplantation if used before the procedure.

CONTRAINDICATIONS

The device must not be used in patients:

- with acute infectious diseases;
- with severe somatic and mental illness;
- with autoimmune and allergic diseases;
- with known allergic or hypersensitive reactions to the components of the product;
- with varices in the area;
- with inflammations or infections in the area of the treatment;
- who tend to develop keloids and hypertrophic scarring;
- who received temporary or permanent implants in the same area;
- who receives anticoagulant treatment;
- who simultaneously receives laser treatment, deep chemical peels or dermabrasion;
- pregnant or breastfeeding women;
- children.

PRECAUTIONS FOR USE

The device must only be used by physicians who have received special training in the mesotherapy injection techniques. The device must be injected into non-inflamed, disinfected, healthy skin. All aseptics and hygiene precautions must be taken to prepare the site, equipment and patient before implantation. Sensitive skin may be pretreated using a local anesthetic patch or cream. After use, discard the syringe, needle and remaining product in a special container.

UNDESIRABLE SIDE EFFECTS

The physician must inform the patient about potential side effects associated with implantation of this device that may occur immediately or within a reasonable delay.

The following events and reactions have been observed with the device or similar products:

- stinging, redness or slight local inflammation which should disappear after 24 to 48 hours;
- mild edema and small hematomas may develop but should disappear within 48 hours.

Patients must inform their physician as soon as possible about any other undesirable effect, or the effects listed above that last longer than the specified duration. The physician should treat these side effects appropriately. Such undesirable side effects must be reported to the distributor and/or manufacturer.

METHOD OF USE

The device must be injected in the dermis of the scalp near hair roots according to the mesotherapy technique.

The device must be implanted aseptically, using single-use sterile syringe which is not supplied by Suissele.

Duration of treatment:

- intensive course: 6 treatments with 2-3 weeks intervals;
- support course: 1 treatment per month.

WARNINGS

- Before use, verify the expiry date and integrity of the package: do not use after the expiry date or if the package is damaged.
- Do not use if the solution is cloudy or if its color differs from that described above.
- Do not mix with other products.
- Do not inject into blood vessels, bones, tendons, ligaments, nerves or muscles.

- Do not overcorrect.
- Once the vial has been opened, use immediately.
- Vial is intended for a single use: do not reuse and do not re-sterilize. Quality and sterility can only be guaranteed for a vial in its original packaging. The re-use of the device creates a potential infection risk for patients or users.

Patients should be advised not to apply any make-up for 12 hours after the injection and to avoid prolonged exposure to the sun, UV rays, frost, as well as using sauna or hammam, for two weeks following implantation.

There are incompatibilities between Sodium Hyaluronate and Quaternary Ammonium Compounds such as Benzalkonium Chloride. Therefore, the device should never be placed into contact with these substances, or with medical-surgical instruments that have been treated with these substances.

STORAGE CONDITIONS

The device should be stored at 5-25°C in a dry place in the original box and protected from direct sunlight, heat or frost. Handle with care.

EXPLANATION OF INTERNATIONAL SYMBOLS	
	Consult instructions for use
	Caution
	Use by
	Batch number
	Do not use if package is damaged
	Sterilized using steam
	Temperature limitation
	Do not re-use
	Keep away from sunlight
	Do not re-sterilize
	Keep dry
	Manufacturer
	Handle with care
	Unique Device Identifier
	Medical Device
	Single sterile barrier system
	CE Mark and Notified Body Number
	Authorized representative in the European Community
	Suissele SA Rue Gallée 6 CH-1400 Yverdon-les-Bains SWITZERLAND Made in Switzerland
	Email: info@suissele.com

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