EN CELLBOCOSTER® SHAPE:
Hyaluronic Acid 18 mg Reshaping Complex

INSTRUCTIONS FOR USE

DESCRIPTION

CELLBOCOSTER® SHAPE: Hyaluronic Acid Reshaping Complex is produced using patented CHAC Technology. It consists of non-crosslinked Hyaluronic Acid with L-Carnitine. The device is designed for reducing local fat accumulation or deposit. It reduces water retention and diminishes puffiness by improving microcirculation and hydration.

The device is a clear, colorless, sterile, injectable solution for single use presented in a case of 6 x 3 mL vials.

COMPOSITION

Each vial contains:

- non-crosslinked Sodium Hyaluronate 18 mg;
- L-Carnitine, Sodium Ascorbyl Phosphate;
- phosphate buffer pH 6.8-7.6 q.s.

INDICATIONS

The device is a resorbable injectable implant for use in the epidermis or dermis for correcting and restoring structure of the skin and subcutaneous tissue to diminish puffiness and cellulite. Improves microcirculation and reduce edema. Helps to prepare for liposuction or surgery and facilitates recovery after them.

CONTRAINDICATIONS

The device must not be used in patients:

- with acute infectious diseases;
- with autoimmune and allergic diseases;
- with known allergic or hypersensitive reactions to the components of the product;
- with varices in the area;
- with inflammations or infections in the area of the treatment;
- who tend to develop keloids and hypertrophic scarring;
- who received temporary or permanent implants in the same area;
- who receives anticoagulant treatment;
- who simultaneously receives laser treatment, deep chemical peels or dermabrasion;
- pregnant or breastfeeding women;
- children.

PRECAUTIONS FOR USE

The device must only be used by physicians who have received special training in the mesotherapy injection techniques. The device must be injected into non-inflamed, disinfected, healthy skin. All asepsis and hygiene precautions must be taken to prepare the site, equipment and patient before implantation. Sensitive skin may be pretreated using a local anesthetic patch or cream. After use, discard the syringe, needle and remaining product in a special container.

UNDESIRABLE SIDE EFFECTS

The physician must inform the patient about potential side effects associated with implantation of this device that may occur immediately or within a reasonable delay.

The following events and reactions have been observed with the device or similar products:

- stinging, redness or slight local inflammation which should disappear after 24 to 48 hours;
- mild edema and small hematomas may develop but should disappear within 48 hours.

Patients must inform their physician as soon as possible about any other undesirable effect, or the effects listed above that last longer than the specified duration. The physician should treat these side effects appropriately. Such undesirable side effects must be reported to the distributor and/or manufacturer.

METHOD OF USE

The device must be injected in the epidermis or dermis of the face, neck, décolleté area, back of the hands, internal face of the arms or body according to the mesotherapy technique.

The device must be implanted aseptically, using single-use sterile syringe which is not supplied by Suisse.

Duration of treatment:

- intensive course: 3 treatments with 2-3 weeks intervals;
- support course: 1 treatment per month.

WARNINGS

- Before use, verify the expiry date and integrity of the package: do not use after the expiry date or if the package is damaged.
- Do not use if the solution is cloudy or if its color differs from that described above.
- Do not mix with other products.
- Do not inject into blood vessels, bones, tendons, ligaments, nerves or muscles.
- Do not overcorrect.
- Once the vial has been opened, use immediately.
- Vial is intended for a single use: do not reuse and do not re-sterilize. Quality and sterility can only be guaranteed for a vial in its original packaging. The re-use of the device creates a potential infection risk for patients or users.
- Patients should be advised not to apply any make-up for 12 hours after the injection and to avoid prolonged exposure to the sun, UV rays, frost, as well as using sauna or hammam, for two weeks following implantation.

- There are incompatibilities between Sodium Hyaluronate and Quaternary Ammonium Compounds such as Benzalkonium Chloride. Therefore, the device should never be placed into contact with these substances, or with medical-surgical instruments that have been treated with these substances.

STORAGE CONDITIONS

The device should be stored at 5-25°C in a dry place in the original box and protected from direct sunlight, heat or frost. Handle with care.

EXPLANATION OF INTERNATIONAL SYMBOLS	
Consult instructions for use	Caution
Use by	Batch number
Do not use if package is damaged	Sterilized using steam
Temperature limitation	Do not re-use
Do not re-sterilize	Keep away from sunlight
Keep dry	Manufacturer
Handle with care	Unique Device Identifier
Medical Device	Single sterile barrier system
Single sterile barrier system with protective packaging outside	CE Mark and Notified Body Number
Authorized representative in the European Community	
Suisse SA Rue Gallié 6 CH-1400 Yverdon-les-Bains SWITZERLAND	Email: info@suisse.com
	Made in Switzerland

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